



## Blood Collection and Processing Instructions for Biobank Samples in the TTM2 and TAME Trials

Version 1.1  
February 26th, 2018



**TABLE OF CONTENTS**

- 1. THE TTM2 AND TAME TRIALS.....3
- 2. BEDSIDE SAMPLE COLLECTION PROCESS.....3
  - 2.1 Blood collection .....3
  - 2.2 Sample collection kits.....3
  - 2.3 Sample collection forms.....4
- 3. SAMPLE PREPARATION AND PROCESSING AT THE SITE LABORATORY .....4
  - 3.1 Materials for sample processing .....4
  - 3.2 How to apply the labels .....4
  - 3.3 Sample processing instruction .....5
  - 3.4 Sample storage log .....6
- 4. SHIPMENT .....6
- 5. CONTACT .....6

APPENDIX A: SAMPLE COLLECTION FORM

APPENDIX B: SAMPLE STORAGE LOG

## 1. THE TTM2 AND TAME TRIALS

The TTM2 (Targeted Temperature Management after Cardiac Arrest 2) and TAME (Targeted Therapeutic Mild Hypercapnia After Resuscitated Cardiac Arrest) trials will collaborate and send samples to a shared biobank located at the Integrated Biobank of Luxembourg (IBBL). Several centres will participate in both trials, but for those patients who only participate in one of the trials the samples will be used in accordance with the consent, by either TTM2 or TAME.

## 2. BEDSIDE SAMPLE COLLECTION PROCESS

### 2.1. Blood Collection

Blood samples will be collected at the site (often the ICU, but could be the Emergency Department or a Ward) from an existing central venous or arterial access device according to the table below.

Blood samples should be drawn at 0h, 24h, 48h, and 72h after randomisation. If laboratory processing according to the instructions is not possible outside “office hours”, samples may be drawn +/- 6 hours from the designated time point. **Efforts should be made to obtain blood samples as close as possible to the ideal collection time points.**

Time point	Blood volume	Collection tube	Processing
0 h*	6 ml 6 ml 2,5 ml	Clot activator tube EDTA-tube PAXgene blood RNA tube	Serum Plasma, Buffy coat -
24 h*	6 ml	Clot activator tube	Serum
48 h*	6 ml 6 ml 2,5 ml	Clot activator tube EDTA-tube PAXgene blood RNA tube	Serum Plasma -
72 h*	6 ml	Clot activator tube	Serum
<b>In total</b>	<b>41 ml</b>		

\*after randomisation

### 2.2. Sample Collection Kits

The biobank IBBL will provide each site sample collection kits for use in patients enrolled in one or both trials. All the materials will be sent to the site study coordinator at the Intensive Care Unit (ICU), who will in turn distribute as appropriate.

The sample collection kits will include pre-labelled collection tubes, Kit ID labels for sample collection forms, and resource materials for the laboratory (cryovials and labels). Each supplied kit contains all the necessary materials for one patient. Each kit will have four bags, one for each time point. After the bedside sample collection, the laboratory materials in the bag should be sent to the laboratory with the samples.

The kit ID, #####-TT-VAR, on the labels, identifies a specific patient and should be recorded in the patient’s medical record. Unused materials must not be used for another trial patient and after the last timepoint the kit should be discarded.

### 2.3. Sample Collection Forms

The TTM2/TAME sample collection form as shown in Appendix A should be used, but may be translated and adapted to local routines. **The site study coordinator should print the sample collection form double sided!** Please follow the instructions in the sample collection form carefully. Complete a new sample collection form for each time point and send it with the samples and the materials for sample processing to the local laboratory. The laboratory will retain each sample collection form and the site study coordinator should collect the forms at regular intervals for data entry purposes.


## 3. SAMPLE PREPARATION AND PROCESSING AT THE SITE LABORATORY

The samples will usually be processed by personnel at the local hospital laboratory, but could be processed by the study site personnel at the ICU if the appropriate training has been given.

### 3.1 Materials for sample processing

The site laboratory will receive blood samples from the ICU four times/patient, or less if a patient is discontinued from the trial. At each timepoint the samples will be sent with a sample collection form and a bag with sufficient materials for the samples to be processed. Each bag will contain cryovials and labels for the cryovials. Boxes for storage of the frozen cryovials will be provided by IBBL before the site commence recruitment.


The ICU will ensure that the correct labels for the cryovials are provided with the blood samples. **The laboratory should be careful to choose the labels for the correct sample type according to the processing instructions.** The labels include the following information (TT being the TTM2/TAME trials identifier):



### 3.2. How to apply the labels

The labels are designed to stick onto clean and dry cryovials at ambient temperature. Use the following method:

- Pull on the support to peel off the label. Be careful not to touch the sticky part of the label.
- Place the long side of the label perpendicular to the cryovial held in the upright direction, as shown in the pictures below. **Both the barcode and the readable text must be visible!**



- Start with placing the white part of the label on the tube, and then rolling it firmly around the cryovial to ensure a good contact and adherence. The transparent part of the label will go over the white part as a second layer.
- Be sure that there are no air bubbles between the label and the cryovial.

If the cryovial for some reason is already frozen, remove the frost on the outer wall before the label is applied. The application procedure must be performed within 4 seconds and the cryovials immediately be put back in the freezer.

### 3.3. Sample processing instructions

Use the cryovials, labels, and sample collection form provided together with the blood collection tubes to process. **The provided labels include the correct kit ID and timepoint, but please be careful to choose the labels for the right sample type (SER, PLA or BFF)!**

#### Serum (at 0h, 24h, 48h, and 72h)

1. Allow the blood in the clot activator tube to clot for 30 min after sampling at room temperature.
2. Label 4 cryovials with the provided labels “#####-TT-VAR-{timepoint}-SER-{01-04}”
3. Centrifuge the clot-activator tube at 2000 g for 10 min at room temperature (slow deceleration at the end of centrifugation).
4. Record the time of centrifugation on the sample collection form.
5. Using a disposable pipette, carefully transfer the serum into an appropriate unused plastic tube. Save the tube containing serum and discard the tube with the clot.
6. Homogenize the serum by up and down movements with a pipette.
7. Aliquot the total volume of serum into the 4 labelled cryovials. Split the serum evenly between the cryovials (if possible at least 500  $\mu$ L/aliquot).
8. **As soon as possible put the labelled cryovials in long term storage at -80°C in the boxes provided.**
9. Record the time of freezing on the sample collection form.

#### Plasma (at 0 h and 48 h)

1. Centrifuge the EDTA- tube as soon as possible at 2000 g for 10 min at room temperature (slow deceleration at the end).
2. Record the time of centrifugation on the sample collection form.
3. Label 4 cryovials with the provided labels “#####-TT-VAR -{timepoint}-PLA-{01-04}”

4. Using a disposable pipette, carefully transfer the plasma into an appropriate unused plastic tube.
5. **Applies only to the sample drawn at 0h (not at 48h):** After removing the plasma, transfer the buffy coat (the thin layer beneath the plasma) and a small portion of the red blood cells with a disposable pipette to a cryovial labelled with "#####-TT-VAR-0-BFF-01". To ensure all the buffy coat is extracted transfer at least 500  $\mu$ L. Discard the collection tube containing the rest of the blood.
6. Homogenize the plasma obtained at point 4 by up and down movements with a pipette.
7. Aliquot the total volume of plasma into the 4 labelled cryovials. Split the plasma evenly between the 4 cryovials (if possible at least 500  $\mu$ L/aliquot).
8. **As soon as possible put cryovials in long term storage at -80°C in the provided boxes**
9. Record the time of freezing on the sample collection form.

#### PAXgene tube (at 0 h and 48 h)

1. Leave the pre-labelled PAXgene tube upright at room temperature for a minimum of 2 hours and a maximum of 24 hours after sample collection.
2. Place the tube in long term storage at -20°C.
3. Record the time of freezing on the sample collection form.

### 3.4 Sample Storage Log

All the TTM2/TAME samples in storage at the local laboratory should be registered in a sample storage log. The cryovials should be placed in the provided storage boxes in the same order as they are listed in the storage log. The log should preferably be completed as soon as the samples are being put in the freezer, but could be created prior to shipment to the biobank IBBL. The sample storage log as shown in Appendix B, or another log with corresponding information, may be used.

The boxes with the cryovials must be stored at -80°C and the PAXgene tubes at -20°C in freezers at the hospital. Temperature logs should be kept for all the freezers storing samples.

## 4. SHIPMENT

Samples should be shipped to IBBL in Luxembourg at the end of the trial, or if necessary at regular intervals (not more than every 6 months to avoid unnecessary shipping costs). The shipment will be organized by the sponsor in collaboration with the site study coordinator. A courier customer code will be given, allowing for the shipment costs to be paid for by the sponsor. Detailed instructions on shipping will be part of a separate form.

## 5. CONTACT

Contact for laboratory/biomarker questions:

TTM2: [helena.levin@med.lu.se](mailto:helena.levin@med.lu.se) or [ttm2@ttm2trial.org](mailto:ttm2@ttm2trial.org)

TAME: [ciara.fahey@ucd.ie](mailto:ciara.fahey@ucd.ie)





## TTM2/TAME SAMPLE COLLECTION FORM

### Sample Processing at the Laboratory

**PLEASE VERIFY THAT THE RED SQUARE ON THE BACK OF THIS FORM HAS BEEN COMPLETED.**  
**If the patient study number or the Kit ID is missing, please contact the ICU!**

Handle the samples according to the provided instructions and complete the information below.

Record date and time in the format DDMMYYYY HH:MM e.g. 15DEC2017 15:30

**Save this form!** The forms will be collected at regular intervals by the site coordinator and the data will be entered in the eCRF at the ICU.

- **Serum** (at 0h, 24h, 48h, and 72h):

**Date and time of centrifugation:**

D	D	M	M	M	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	H	H	:	M	M

**Date and time placed in freezer:**

D	D	M	M	M	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	H	H	:	M	M	

- **Plasma** (at 0 h and 48 h):

**Date and time of centrifugation:**

D	D	M	M	M	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	H	H	:	M	M	

**Date and time placed in freezer:**

D	D	M	M	M	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	H	H	:	M	M	

- **PAXgene tube** (at 0 h and 48 h):

**Date and time placed in freezer:**

D	D	M	M	M	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	H	H	:	M	M	

**Samples handled by** (name/initials):

.....

If any deviations occurred, please specify:





## TTM2/TAME SAMPLE STORAGE LOG

Patient's Kit ID #####-TT-VAR (first 8 digits on the label)	Time Point	Sample Type	Number of Vials in Storage	Storage in Box No(s)	Date and Initials
	0	Serum			
	0	Plasma			
	0	Buffy coat			
	0	PAXgene			
	24	Serum			
	48	Serum			
	48	Plasma			
	48	PAXgene			
	72	Serum			
	0	Serum			
	0	Plasma			
	0	Buffy coat			
	0	PAXgene			
	24	Serum			
	48	Serum			
	48	Plasma			
	48	PAXgene			
	72	Serum			
	0	Serum			
	0	Plasma			
	0	Buffy coat			
	0	PAXgene			
	24	Serum			
	48	Serum			
	48	Plasma			
	48	PAXgene			
	72	Serum			